

“witeg” **Wisetherm**[®] **WHM** Aluminum-case Beaker Heating Mantles with Controller

Nickel Chrome Heating Element, Built-in Temperature Controller, 100 ~ 5.000 ml, 450 °C

- Full Choices of Vessel, Capacity, Housing and Control are Offered!
- **Nickel Chrome Heating Element - Coil Type**
: Acid-Resistant, Long Life Time
- Built-in Temperature Controller
- All Kinds of Round Flasks Can be used - DURAN/Kimax/Pyrex/Iwaki
- CE Certified

Catalog No.		Flask Capacity (ml)	Flask Dia. ϕ (mm)	Watts (W)	Inside Depth (mm)	Mantle Outside $\phi \times H$ (mm)	Support Rod $\phi \times H$ (mm)	Application (Flask-brand)
with 230 V	with 110 V							
DH.WHM12141	DH.WHM12151	100	50 ~ 52	80	40	ϕ 145 \times H 105	ϕ 13 \times H 440	Common-use
DH.WHM12142	DH.WHM12152	250	68 ~ 70	140	65	ϕ 145 \times H 140	ϕ 13 \times H 440	Common-use
DH.WHM12143	DH.WHM12153	400	79 ~ 81	240	70	ϕ 145 \times H 140	ϕ 13 \times H 440	for DU.
DH.WHM121431	DH.WHM121531	400	76 ~ 78	240	75	ϕ 145 \times H 140	ϕ 13 \times H 440	for Iw/Ki/Py.
DH.WHM12144	DH.WHM12154	600	88 ~ 90	240	75	ϕ 145 \times H 140	ϕ 13 \times H 440	Common-use
DH.WHM12145	DH.WHM12155	800	98 ~ 100	300	90	ϕ 190 \times H 165	ϕ 13 \times H 440	Common-use
DH.WHM12146	DH.WHM12156	1.000	105 ~ 108	350	90	ϕ 190 \times H 165	ϕ 13 \times H 440	Common-use
DH.WHM12147	DH.WHM12157	2.000	130 ~ 132	450	125	ϕ 270 \times H 200	ϕ 13 \times H 660	Common-use
DH.WHM12148	DH.WHM12158	3.000	151 ~ 153	550	160	ϕ 270 \times H 235	ϕ 13 \times H 660	for DU/Iw.
DH.WHM121481	DH.WHM121581	4.000	160 ~ 162	700	160	ϕ 270 \times H 235	ϕ 13 \times H 660	for Ki/Py.
DH.WHM12149	DH.WHM12159	5.000	169 ~ 171	800	160	ϕ 270 \times H 235	ϕ 13 \times H 660	for DU.

※ Remarks : DU. : DURAN Ki. : Kimax Iw. : Iwaki Py. : Pyrex



“witeg” **Wisetherm**[®] **WHM** Aluminum-case Stirrer Beaker Mantles with Controller

Nickel Chrome Heating Element, Stirring & Temperature Controller, 100 ~ 5.000 ml, 450 °C

- Full Choices of Vessel, Capacity, Housing and Control are Offered !
- **Nickel Chrome Heating Element - Coil Type**
: Acid-Resistant, Long Life Time
- Built-in Stirring & Temperature Controller
- All Kinds of Round Flasks Can be used - DURAN/Kimax/Pyrex/Iwaki
- CE Certified

Catalog No.		Flask Capacity (ml)	Flask Dia. ϕ (mm)	Watts (W)	Inside Depth (mm)	Mantle Outside $\phi \times H$ (mm)	Support Rod $\phi \times H$ (mm)	Application (Flask-brand)
with 230 V	with 110 V							
DH.WHM121601	DH.WHM121611	100	50 ~ 52	80	45	ϕ 145 \times H140	ϕ 13 \times H440	Common-use
DH.WHM121602	DH.WHM121612	250	68 ~ 70	140	70	ϕ 190 \times H165	ϕ 13 \times H440	Common-use
DH.WHM121603	DH.WHM121613	400	79 ~ 81	240	70	ϕ 190 \times H165	ϕ 13 \times H440	for DU.
DH.WHM121604	DH.WHM121614	400	76 ~ 78	240	70	ϕ 190 \times H165	ϕ 13 \times H440	for Iw/Ki/Py.
DH.WHM121605	DH.WHM121615	600	88 ~ 90	240	72	ϕ 190 \times H165	ϕ 13 \times H440	Common-use
DH.WHM121606	DH.WHM121616	1.000	105 ~ 108	350	105	ϕ 270 \times H200	ϕ 13 \times H440	Common-use
DH.WHM121607	DH.WHM121617	2.000	130 ~ 132	450	140	ϕ 270 \times H235	ϕ 13 \times H660	Common-use
DH.WHM121608	DH.WHM121618	3.000	151 ~ 153	550	140	ϕ 270 \times H235	ϕ 13 \times H660	for DU/Iw.
DH.WHM121609	DH.WHM121619	4.000	160 ~ 162	700	170	ϕ 325 \times H270	ϕ 13 \times H660	for Ki/Py.
DH.WHM121610	DH.WHM121620	5.000	169 ~ 171	800	170	ϕ 325 \times H270	ϕ 13 \times H660	for DU.

※ 2 pcs of PTFE Magnetic Stirrer Bar Offered : Oval-type, ϕ 6,0 \times L 15 mm (CW.001.615) and ϕ 10 \times L25 mm (CW.001.625.10) one each

※ Remarks : DU. : DURAN Ki. : Kimax Iw. : Iwaki Py. : Pyrex