

Name	Value
Pharmacy refrigerator	MD40DD-P600SEQMXXXXX1SN4B
EAN code	4770104040054
Gross capacity, l	385
Net capacity, l	270
Factory temperature setting, °C	+5
Fluctuation, °C	0,5
Gradient, °C	2,9
Increase in temperature from +5 °C to +10 °C (empty, +25 °C TU)	40 min.
Decrease in temperature from +25 °C to +5 °C (empty, +25 °C TU)	40 min.
Ambient temperature range, °C	from +10 to +35
Airborne acoustical noise emissions, dB(A) re1 pW	50
Refrigerant	R600a
Daily consumption, kWh/24h	1,440
Annual energy consumption, kWh/year	525
Rated voltage, V	220 - 240
Rated frequency, Hz	50
Self closing door	Yes
Access port	Yes
Cooling system	Fan forced cooling
Defrost	Automatic
Volt-free alarm contact	Yes
Interior light	LED ceiling lighting
Malfunction warning signal	Optical and acoustical
Power failure alarm	Immediately upon power failure for min. 12h
Connection to the external alarm system	Yes
Data logger	Integrated, USB interface
Protection against low temperature	Yes
Number of shelves	5 adjustable +1 bottom
Shelf size (W x D), mm	380 x 495, bottom 445 x 180
Maximum load per shelf, kg	40
Shelf material	Wire shelves coated by white polymer paint
Side wall material	Steel, white painted
Inner cabinet material	HIPS polystyrene
Type of door	Full-panel door
Type of lock	Mechanical
Adjustable foot at the front	2
Wheels behind	2

Name	Value
External dimensions (H x W x D), mm	2025 x 600 x 620
Internal dimensions (H x W x D), mm	1645 x 500 x 418
Package dimensions (H x W x D), mm	2080 x 620 x 630
Weight net, kg	80.5
Weight gross, kg	86