

Data sheet

Dygestorium Compact Line DCL 15.00 SS WALK-IN



Standard equipment:

- illumination of working chamber
- air-flow sensor (Schneider FM 550-A-0-E)
- Declaration of conformity CE: PN-EN 14175; PN-EN 14727; PN-EN 14056; PN-EN 13150

Additional options:


- electrical sockets
- cold water tap

TECHNICAL DATA - DIMENSIONS
overall dims [mm]

width A	1500
height B	2385-2870
depth C	1200

working space dims (mm)

width D	1170
window width D1	1110
height E	2145
depth F	845


TECHNICAL DATA

recommended airflow [m ³ /h]	750 - 1200
nominal power [W]	82
voltage**	230 VAC 50Hz
diameter of the air outlet stub [mm]	200
airflow sensor (Schneider FM 550-A-0-E)	compliance with PN-EN 14175-2
electrical insulation	class I
required capacity of the ventilation system	0,3 - 0,5 [m/s]
standard ventilation system	double rear wall
frame and housing material	The frame is made entirely of galvanized steel. External housing made of galvanized and painted with epoxy paint steel
working chamber wall finishing	galvanized steel painted with epoxy paint