

Name	Value
Pharmacy refrigerator	MD40DD-P301BEQMXXXXX1SN4B
EAN code	4770104040115
Gross capacity, l	385
Net capacity, l	290
Factory temperature setting, °C	+5
Fluctuation, °C	0,5
Gradient, °C	2,9
Increase in temperature from +5 °C to +10 °C (empty, +25 °C TU)	40 min.
Decrease in temperature from +25 °C to +5 °C (empty, +25 °C TU)	40 min.
Ambient temperature range, °C	from +10 to +35
Airborne acoustical noise emissions, dB(A) re1 pW	44
Refrigerant	R600a
Daily consumption, kWh/24h	1,280
Annual energy consumption, kWh/year	467
Rated voltage, V	220 - 240
Rated frequency, Hz	50
Self closing door	Yes
Access port	Yes
Cooling system	Fan forced cooling
Defrost	Automatic
Volt-free alarm contact	Yes
Interior light	Vertical LED light strip on right side
Malfunction warning signal	Optical and acoustical
Power failure alarm	Immediately upon power failure for min. 24h
Connection to the external alarm system	Yes
Alarm indicator on screen	Visible until removed manually
Data logger	Mobile app connected via bluetooth
Protection against low temperature	Yes
Number of shelves	5 adjustable + 1 bottom
Shelf size (W x D), mm	380 x 495, bottom 445 x 180
Maximum load per shelf, kg	40
Shelf material	Wire shelves coated by white polymer paint
Side wall material	Steel, white painted
Inner cabinet material	HIPS polystyrene
Type of door	Double panel selective glass in aluminium profile
Type of lock	Mechanical
Adjustable foot at the front	2

Name	Value
Wheels behind	2
External dimensions (H x W x D), mm	2025 x 600 x 620
Internal dimensions (H x W x D), mm	1645 x 500 x 418
Package dimensions (H x W x D), mm	2080 x 620 x 630
Weight net, kg	80,5
Weight gross, kg	86